

meywalk[®] 2000

The springy walker



MEYWALK[®] 2000 is supplied standard with seat trunk support, rear stop, handlebar, fender wheels and parking brakes.

Meywalk[®] 2000 Freedom of movement for children

MEYWALK[®] 2000 from Meyland-Smith is the improved version of the famous Meywalk[®] walking support. Now it has been stepless quick adjustment of body support and seat height.

MEYWALK[®] 2000 is available in three sizes. For children from approx. 6 years and up to adults.

Seat and trunk support has spring suspension helpful for a better and more natural time settlement. The seat unit can easily be lowered, making entry and exit easier for the user and less stressful of helpers.

The user's posture continuously adjustable from vertical to 15 ° stooping, and that makes it possible to find the body inclination that gives the user best control.

Meywalk[®] 2000 is made of best quality materials ensures long life, maximum reliability and minimal maintenance.

Meywalk[®] 2000 is therefore suitable for both indoor and outdoor use.



Dejay Medical & Scientific Pty Ltd

Telephone: (02) 9838 8869

www.dejay.com.au

meywalk[®] 2000

Accessories

- Anti-tip supports
- Extra-long spacer tubes (Provides 6 inches in height - only size Large)
- Height Decreases bracket (Gives 8 cm less in height)
- Drag Brakes
- Non-reverse brakes (can be deactivated)
- Leg guides
- Swivel lock system for front wheels (selectable)
- Hip pads
- Crossbar – no seat option
- Leg separation plate



Meywalk 2000 Large

MEYWALK[®] 2000's seat and trunk support has spring suspension, and springs this is available in four different strengths, which should fit the user's weight:

- ∴ Min. 20 kg = White springs.
- ∴ Min. 45 kg = Yellow springs.
- ∴ Min. 65 kg = Blue springs.
- ∴ Min. 80 kg = Red springs.

Body Support size must be disclosed when ordering!

MEYWALK[®] 2000 has a built in raising-lowering device facilitating. This makes it easier for the user to get in or out of the walking support. The lock buttons in front of Meywalk[®] 2000 are released and the swing bar tilted up, whereby the seat unit is lowered approx. 15 cm in size Medium and approx. 25 cm for size Large. When the user is in place, raise seat unit again by pressing the swing bar down until it locks with an audible click.



Meywalk 2000 Medium



Meywalk 2000 Small

Technical data:	Small	Medium	Large
Length:	85 cm	95 cm	95 cm
Width:	68 cm	70 cm	70 cm
Seat/groin height:	48 - 66 cm	62 - 86 cm	76 - 100 cm
Trunk Support Height:	81 - 100 cm	94 - 117 cm	108 - 131 cm
Wheel size, front:	150 x 50 mm	200 x 50 mm	200 x 50 mm
Wheel size, rear:	400 x 50 mm	400 x 50 mm	400 x 50 mm
Weight:	22 kg	25 kg	26 kg
Max. user weight:	100 kg	100 kg	100 kg



Dejay Medical & Scientific Pty Ltd

Telephone: (02) 9838 8869

www.dejay.com.au