



Designed and hand built in Germany.





Walking aid FLUX, size 2

### Pictured configuration:

Durable, compact frame in a light aluminum construction.

Grab rail height- and width adjustable. Black and flat-free, polyurethan front swivel wheels. Black and flat-free, polyurethan rear wheels. Engageable and disengageable free-wheeling system on rear wheels. Anit-tips without wheels. Frame colour: luminous yellow.



### Natural walking.

### Natural walking with the *FLUX* walking aid:

**FLUX** is a walking aid that promotes natural walking. Opening out in walking direction and optically and physically free from interfering structural parts **FLUX** encourages upright posture and natural walking movements.







FLUX size 2.



Front view.



Front view.

Extremely compact folding size.

FLUX size 4.

New frame colours. All options available online.

**FLUX** has an extremely compact folding size and is available in four sizes with a selection of trendy colour options. All **FLUX** accessories including pictures can be found online: **www.hoggi.de** 





Walking aid **FLUX**, size 3

### Pictured configuration:

Durable, compact frame in a light aluminum construction. Grab rail with universal grips. Black and flat-free, polyurethan front swivel wheels. Black and flat-free, polyurethan back wheels. Engageable and disengageable free-wheeling system. Anti-tips without wheels. Frame colour: luminous red.



# The perfect position - in every situation.

### FLUX optimally adapts to your body height: Adjustable grip bar variations.

Height-adjustable grip bars make it easy to adjust to different body heights. The grip bar with universal grips (optional) allows an adjustment of the grip width.







Vertical adjustment.



Horizontal adjustment.



Adjustable grip width.

## Stepless height adjustment. Safety aspect: Free-wheeling system.



Free-wheeling system sizes 1 & 2.



Free-wheeling system size 3.



Free-wheeling system size 4.

### Safety aspect: Free-wheeling system.

Innovative and reliable reverse-roll locks ensure your safety and prevent unintended backward-rolling of your *FLUX*.

# FLIJX

# A selection of options.



Anti-tips without wheels.



Anti-tips with wheels.



Dynamic hip control.



Swivel lock.



Seat folded down.



Seat folded up.

Safety aspect: Anti-tips. Dynamic hip control and resting seat.

The sturdy anti-tips keep you safe and prevent your *FLUX* from tipping over.

Anti-tips with or without wheels is available for all sizes.

The hip control is multidirectionally adjustable and the hip belt offers additional support while walking. The seat can be easily folded up or down according to your needs.



Technical data	Size 1	Size 2	Size 3	Size 4
Grab rail standard				
Grip height	48 - 61 cm	57 - 81 cm	69 - 87 cm	82 - 102 cm
Grip width	41 cm	46 cm	51 cm	51 cm
Grab rail with universal grips				
Grip height	43 - 68 cm	51 - 88 cm	62 - 94 cm	75 - 102 cm
Grip width	29 - 41 cm	34 - 46 cm	39 - 51 cm	39 - 51 cm
Lower arm support				
Height	48 - 69 cm	57 - 89 cm	69 - 95 cm	82 - 110 cm
Distance forearms	29 - 41 cm	34 - 46 cm	39 - 51 cm	39 - 51 cm
Seat height	37 cm	44 cm	54 cm	58 cm
Frame				
Length	65 cm	80 cm	88 cm	90 cm
Width	64 cm	69 cm	74 cm	78 cm
Wheel diameter				
Front	15 cm	15 cm	15 cm	20 cm
Back	20 cm	20 cm	20 cm	30 cm
Load capacity	60 kg	75 kg	75 kg	100 kg
Weight	6,5 kg	6,8 kg	7,1 kg	9 kg



# SCAN ME to view the product video in HD quality!

Individual custom-made solutions.

Customize your *FLUX* according to your own medical and physical needs. Everything is possible. The custom-made solutions are designed especially to your individual needs and guarantee perfect patient-centered care. Let us know - we are happy to help!





Dejay Medical 9 Anvil Road Seven Hills NSW 2147 Australia

P: 02 9838 8869 E: info@dejay.com.au W: www.dejay.com.au

#### HOGGI

Partnering together to mobilse kids.

### **HOGGI** GmbH

- Eulerstraße 27D-56235 Ransbach-Baumbach
- Telefon: (+49) 2623 / 92 499-0 Telefax: (+49) 2623 / 92 499-99

www.hoggi.de